

ENT Single Use Suction Range



Tips



Handles



Speculum

Welcome to GBUK Healthcare

ENT single use suction specialists

METAL SUCTION TIPS

- Single use, for surgery and examination of the ear
- ETO sterilised, packed in a peel pack
- Curved angle tip specifically designed for microscopic examination and surgery
- Medical grade stainless steel, luer slip handle




“ A comprehensive range of metal suction tips to suit all ENT procedures.”

Magnified showing the comparative quality of metal suction tip ends

NEW METAL SUCTION TIP STAPES

- Tube in tube for rigidity and fine tip end lumen
- Single use, sterile, packed in a peel pack
- Special suction tip for stapes surgery



SUCTION TIP HANDLES


- Single use
- With or without suction control, should be used as a connector between suction tip and suction tubing
- Sterile, single packed in a peel pack
- White ABS plastic



“Handles that allow more controlled suctioning.”

TRANSPARENT HANDLE SUCTION TIPS (Frazier/Ferg)


- Single use, sterile, packed in a peel pack
- For surgical operations
- Medical grade stainless steel, transparent handle (PMMA)



“A comprehensive range of products for all of your ENT suction needs.”

MYRINGOTOMY SCALPEL


- Single use, sterile, single packed in a peel pack
- Triple cut edge made of medical grade stainless steel
- Used for paracentesis



“An ENT range with a clear understanding of procedural needs.”

MYRINGOTOMY KIT

- Single use
- Readymade kit for paracentesis surgery
- Sterile, single packed in a peel pack



“Convenient, fast and safe.”

ANCILLARY

“High quality accessories complementary to the ENT suction range.”

MEDIBENDER



- A bending instrument for mediplast suction tips
- Stainless steel, non sterile
- Can be autoclaved

EAR FLUSHING TIP



- Non-sterile, single packed in a peel pack
- Fits on catheter tip syringe
- Safety wing to prevent ear trauma

NASAL SPECULUM



- Single use, single packed in a peel pack
- Reversed spring effect allows physician to use both hands.

AURAL SPECULUM



- Oval design of cone and side cut end for improved sight
- For examination procedures, fits Peter's and Siegler's pneumatic glass

ANTRUM IRRIGATION SET



- Single use
- 1,4 x80 mm needle, 200mm extension tube with luer lock connector
- Sterile, single packed in a peel pack



www.gbukhealthcare.com

Other GBUK Healthcare surgical products



MGE



Serres



TLC Retractor

+44 (0)1757 288587

Blackwood Hall Business Park, North Duffield, Selby, North Yorkshire, YO8 5DD
Tel: +44 (0)1757 288587 Fax: +44 (0)845 280 3160 Email: info@gbukhealthcare.com

www.gbukhealthcare.com